



Gel "E" 2" General Use Gel/Foam Wheelchair Cushion

Gel-U-Seat™ Lite General Use 2" Gel/Foam Wheelchair Cushion

Item #	Product Dimensions	Product HCPCS	UOM
14888	18" (W) X 16" (D) X 2" (H)	E2601	1/ea

- The Medicare approved General Use 2" gel-foam wheelchair cushion will provide comfort and support for the wheelchair bound patient.
- High-Density foam encapsulates a viscous gel bladder that provides proper immersion and envelopment for pressure redistribution.
- The cover is fluid resistant low-shear nylon; zippered and removable with a durable vinyl bottom.
- Carry handle included for convenient transport.
- Meets flammability standard TB117.

specifications

Cover Material	Fluid-Resistant Nylon Waterproof Vinyl Bottom
Weight Capacity	250 lbs.
Warranty	1 Year Limited

Item #	Product Dimensions	Product HCPCS	UOM
8040-1	16" (W) X 16" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-2	18" (W) X 16" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-5	18" (W) X 18" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-3	20" (W) X 16" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-4	20" (W) X 18" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-10	20" (W) X 20" (D) X 2" (H)	E2601	1/ea, 4/cs
8040-6	22" (W) X 16" (D) X 2" (H)	E2602	1/ea, 4/cs
8040-7	22" (W) X 18" (D) X 2" (H)	E2602	1/ea, 4/cs
8040-9	24" (W) X 18" (D) X 2" (H)	E2602	1/ea, 4/cs

- Offering great versatility, the Gel-U-Seat™ Lite provides comfort to almost all wheelchair users.
- Unique dual-chambered gel bladder keeps the pelvis in contact with pressure redistributing cooling gel.
- The cover is fluid resistant low-shear nylon; zippered and removable with a durable vinyl bottom.
- Secures easily to wheelchair with 12" tie straps.
- Meets flammability standard TB117.

specifications

Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
Weight Capacity	Up to 20" 275 lbs., 22" and up 500 lbs
Warranty	1 Year Limited