

Group II Powered Overlay



14427 - Quick Inflate Pump  
(Sold Separately)

## Med-Aire 5" Alternating Pressure and Low Air Loss Overlay System

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
14025N	Alternating Pressure/Low Air Loss Overlay System	36" (W) X 80" (L) X 5" (H)	E0372	1/ea

- This combination therapy overlay system provides both alternating pressure and low air loss to optimize pressure redistribution and manage skin maceration.
- Designed to prevent, treat and manage pressure ulcers in the home, long term care, or acute care setting.
- Pump delivers air volume of 8 LPM. (Figure A)
- The nylon low shear cover material is water resistant, vapor permeable, machine washable and secures to the overlay with straps.
- Each one of the 17, 5" deep air bladders is easily removed and replaced.
- Pillow function maintains the 3 air cells at the head of the bed in static mode for patient comfort.
- CPR valve allows for rapid deflation.
- EZ Lock/Quick Release Tubing Connector has pressure tabs that allow easy and quick disconnect. (Figure B)
- Adjustable patient comfort settings.

### specifications

Control Unit Airflow	8 LPM
Control Unit Cycle Time	10 Minutes
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant Nylon
Weight Capacity	350 lbs.
Warranty	1 Year Limited