



Customize Advanced Digital Pump to Individual Patient's Weight



Standardize Inventory: One Pump for 3 Different Support Surfaces



14427 - Electric Pump (sold Separately)  
Replaces Mason # 8151-EP

**drive Mason™** Masonair AS8800 8" Alternating Pressure and Low Air Loss Mattress System

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
AS8800	Alternating Pressure/Low Air Loss Mattress System	80" (L) x 36" (W) x 8" (H)	E0277	1/ea
AS8800-84	84" Alternating Pressure Power Support Surface Mattress	84" (L) x 36" (W) x 8" (H)	E0277	1/ea
AS8800-RR	80" Raised Rail Alternating Pressure Power Support Surface Mattress	80" (L) x 36" (W) x 10" (H)	E0277	1/ea

- The microprocessor sensor technology provides continuous feedback to equalize pressure and maximize patient comfort.
- The mattress consists of (20) 8" high cell-on-cell designed air bladders.
- In case of a power failure the bottom 4" air cell will remain inflated offering sufficient weight support for the patient.

- The 10 adjustable comfort settings allows pressure to be customized by patient weight.
- Static Mode can suspend alternating mode and has auto-recovery after 30 minutes.
- Pillow Top feature keeps head and neck stable.
- Auto-Firm feature for max inflation when transferring patient and has auto-recovery after 30 minutes.
- Seat Inflation mode offers additional support when head of bed is raised for an extended period of time.
- Visible and audible alarms with Alarm Mute feature (3 minute auto-return feature).
- CPR valve allows for rapid deflation.

**specifications**

Pump Airflow	8 LPM
Pump Cycle Time	10 , 15, 20, 25 Minutes
Pump Type	Compressor
Pump Alarms	Visual / Audible
Mattress Dimensions	(80", 84") x 36" x 8"
Mattress Cover Material	Fluid-Resistant Nylon
Weight Capacity	350 lbs.
Warranty	1 Year Limited